


Finnigan™ ELEMENT2



**High Performance
High Resolution ICP-MS**

Finnigan™ ELEMENT2

High Performance High Resolution ICP-MS

An Instrument for every Application

A success story around the world with over 300 instruments installed, for example:



Belgium: OCAS.

Research centre for the application of steel with New Wave laser ablation system.

www.ocas.be



Sweden: Analytica AB.

Scandinavia's leading contract lab with Emma Engström.

www.analytica.se



Canada: Seastar Chemicals Inc.

A manufacturer of high purity chemicals with Dr. Brad McKelvey.

www.seastarchemicals.com



USA: Micron Technology Inc.

Advanced semiconductor solutions with Kevin Coyle.

www.micron.com



Germany: Institute for Transuranium Elements.

Shown with a glovebox for radioactive sample measurements.

<http://itu.jrc.cec.eu.int>



USA: Desert Research Institute.

Hydrologic Sciences Division with Steve Lambert.

www.dri.edu

- Multielement analysis across the periodic table covering a mg/L to sub pg/L concentration range
 - *Compatible with inorganic and organic solution matrices and solids*
- High mass resolution to access spectrally interfered isotopes
 - *Produces unambiguous elemental spectra*
- A multielemental detector for transient signals
 - *For example, CE, HPLC, GC, FFF and laser ablation*
- High precision isotope ratios
 - *Independent of interferences or interfered isotopes*
- Fully automated tuning and analysis
 - *In conjunction with a comprehensive, customizable quality control system*
- Reliability and robustness to serve as a 24 / 7 production control tool
 - *Highest sample throughput*
- Highest flexibility and accessibility to serve as an advanced research tool




High Mass Resolution

The unequivocal separation of analyte ions from spectral interferences is a prerequisite of accurate and precise analysis. High mass resolution is the universal means for this separation.


Spectral interferences are the main limitation of ICP-MS. The argon plasma gas, water, acid and the sample matrix itself can combine to introduce a wide range of polyatomic ion species. The resultant interfering species may have the same nominal mass as an analyte ion and thus return a falsely high value for the analyte. Numerous strategies have been used in order to try to minimize or circumvent the formation of these spectral interferences. These include mathematical corrections, special sample introduction systems, special plasma parameters and collision/dynamic reaction cells to neutralize part of the interferences. High resolution with a sector-field mass spectrometer simply distinguishes the analyte from interference by difference in mass.

The capability of high mass resolution is a feature unique to the Finnigan™ ELEMENT2. This ability can be used for quantification and isotope ratio analysis for nearly the whole periodic table, and in almost all matrices.


Even in a sample matrix as simple as ultra pure water (UPW), interferences do exist especially at low analyte concentration levels:




Phosphorus in UPW, Medium Resolution



Silicon in UPW, Medium Resolution



Calcium in UPW, Medium Resolution




Zinc in H_2SO_4 (10 % w/w), High Resolution

The more complex the sample matrix, the wider the range of interferences that will occur. The main advantage of high mass resolution as a technique to remove spectral interferences is that it is not just limited to a particular type of interference. With high mass resolution, iron at m/z 56 is easily separated in a simple matrix from $^{40}\text{Ar}^{16}\text{O}$ as well as in more complex matrices from, for example, $^{40}\text{Ca}^{16}\text{O}$.


Arsenic at m/z 75 can be determined in HCl, separated from the $^{40}\text{Ar}^{35}\text{Cl}$ interference, or, in a matrix containing Ca and chloride, from $^{40}\text{Ca}^{35}\text{Cl}$.

Matrices such as mineral acids and organic solvents, which form a plethora of spectral interferences, can be easily analyzed using the high resolution mode.

High resolution results in simple clear spectra and does not create new interferences.



Iron in Liquid Crystal Matrix, Cold Plasma, Medium Resolution



Nickel Isotopes in Groundwater, Medium Resolution

Resolution Specification

(10 % peak valley definition)

3 fixed resolutions:

Low Resolution > 300

Medium Resolution > 4000

High Resolution > 10000

Principle

The Finnigan ELEMENT2 is a double focusing magnetic sector field ICP-MS.

Plasma and Interface

The argon plasma ion source and sampling interface of the Finnigan ELEMENT2 are at ground potential¹. This enables the straightforward coupling of peripherals like HPLC, CE, GC and laser ablation. The Finnigan ELEMENT2 interface reduces the initial kinetic energy spread from ~20 to ~5 eV by capacitively decoupling the plasma from the load coil, using a grounded guard electrode. This reduced energy spread increases the ion transmission and delivers superior sensitivities at all resolutions.

The ion transfer optics focus the ions from the plasma interface on the entrance slit of the double focusing analyzer. The Finnigan ELEMENT2 ion transfer optics are designed for low background, highest sensitivity and minimum mass bias at maximum stability².

High Resolution

The Finnigan ELEMENT2 is able to fully automatically change between three fixed resolutions by switching the positions of the entrance and exit slits in < 1 s. The patented design of the fixed slit mechanism³ offers maximum stability and reproducibility of resolutions.

¹ Patent issued: US 5552599, GB 2282479

² Patent issued: US 5625185

³ Patent issued: US 5451780, GB 2281438



- Highly laminated, water cooled magnet, m/z : 7 - 240 - 7 in < 150 ms
- 25 ppm/8 hour mass stability

MR LR HR

- 3 fixed resolutions:
 R = 300
 R = 4000
 R = 10000




- Ion Transfer Optics
 Guarantees a flat response curve,
 low background and high sensitivity



- Simultaneous measurement in analog and counting modes
- $> 10^9$ linear dynamic range
- < 0.2 cps dark noise
- Fully automatic cross calibration



- Interface at ground potential
- Easily changeable cones
- Resistant against chemicals



- Torch with Guard Electrode (GE)
The GE decreases the ion energy spread, thus increasing ion transmission. This and the high acceleration voltage used in a sector field ICP-MS results in increased sensitivity. The GE is also required for Cold Plasma measurements

Mass Separation

The magnetic field disperses ions according to their mass and energy. The magnet used in the Finnigan ELEMENT2 is specifically designed for use in ICP-MS applications. It is relatively small (sufficient for the mass range 0 – 260 u), highly laminated and efficiently water-cooled for the highest mass stability. Changing of the magnetic field is controlled by a magnetic field regulator with a new high power stage, which delivers the fastest scan speed ever possible with a magnetic sector field instrument⁴.

After passing through the magnetic field the ions enter the electrostatic analyzer for energy focusing. The combination of the magnetic and electrostatic fields results in the double focusing, high resolution properties of the Finnigan ELEMENT2.

Detection System

The Finnigan ELEMENT2 is equipped with a discrete dynode detector system. Rather than having the ion beam directly strike the detector to initiate an electron cascade, the secondary electron multiplier implemented in the Finnigan ELEMENT2 uses a conversion dynode at -8 kV, producing a uniform response across the mass range. The detector is linear over nine orders of magnitude – from ppq to ppm concentrations. The quantification of trace and major elements is therefore possible in a single analysis.

⁴ Thermo Electron Application Note AN30011_E12/03C

Resolution

Low Resolution


Low resolution ($R = 300$) is used for the analysis of non-interfered isotopes. In this mode the Finnigan ELEMENT2's sensitivity is the highest of all commercially available ICP-MS instruments. Additionally, the flat top peak shape is an advantage for high precision isotope ratio measurements.

Medium Resolution


Medium resolution ($R = 4000$) guarantees interference-free analysis for most elements in the majority of sample matrices. For example, transition elements are routinely measured in medium resolution due to the formation of many interfering polyatomic species in the mass range 24 – 70 u.

High Resolution


High resolution ($R = 10000$) is used for the analysis of elements in the most challenging sample matrices. For example, high resolution is used to separate As and Se from argon dimer interferences and argon chloride interferences in chlorine matrices, heavy rare earth elements from light rare earth element oxides in geological matrices, platinum group elements from argon-transition metal molecular species, and/or the oxides of Hf, Ta and W.




100 ng/L ^{115}In , Low Resolution



100 ng/L ^{115}In , Medium Resolution



100 ng/L ^{115}In , High Resolution




100 µg/L Hafnium, 100 ng/L Platinum, High Resolution


Since the change in mass resolution is achieved by changing the width of the entrance and exit slits of the mass spectrometer, the instrumental sensitivity of a high resolution ICP-MS is dependent on the resolution mode used. Therefore, the Finnigan ELEMENT2 with three fixed resolutions has three sensitivities: the wider the slit, the higher the sensitivity.

Fixed sensitivity ratio between resolutions: independent of mass and matrix.

Even in high resolution mode, the intrinsic sensitivity of the Finnigan ELEMENT2 provides sub ppt detection limits.



Chromium in Blood, Medium Resolution



Thulium in Europium Matrix, High Resolution

Sensitivity Specification

Low Resolution (R = 300)

$^{115}\text{In} > 1 \times 10^6$ cps/ppb

Medium Resolution (R = 4000)

$^{115}\text{In} > 1 \times 10^5$ cps/ppb

High Resolution (R = 10000)

$^{115}\text{In} > 1.5 \times 10^4$ cps/ppb

Background in all 3 resolutions < 0.2 cps

Sensitivity and Stability

Sensitivity and Low Background


It is obvious that the highest instrumental sensitivity is essential to achieve the lowest detection limits. However, it is the signal to noise ratio that dictates the detection limit. The Finnigan ELEMENT2 guarantees an off-peak background of < 0.2 cps for all three resolutions. Detection limits in the fg/L range are possible.

Detection Limits in Complex Matrices

Even in cases where the lowest detection limits are not the goal, the exceptionally high sensitivity of the Finnigan ELEMENT2 offers an important advantage. With the higher sensitivity, higher dilution factors for complex matrices can be used without sacrificing detection limits.

Sensitivity and 'Dilute and Shoot'


The requirement for matrix separation for samples such as seawater, mineral acids or organic solvents is replaced by a simple dilution. Using this 'dilute and shoot' approach, the matrix load on the sample introduction system, plasma and interface is reduced. This is also particularly important for the analysis of nuclear sample matrices where the minimization of waste is of paramount importance.




Mineral Water, undiluted (5 % HNO₃), ~200 fg/L Radium, Low Resolution

	<i>Resolution</i>	<i>LoD [ng/L] in solution</i>
²³ Na	Low	1.3
⁹⁰ Zr	Low	0.06
¹⁰⁷ Ag	Low	0.03
¹⁵¹ Eu	Low	0.04
⁴⁷ Ti	Medium	0.3
⁵² Cr	Medium	0.3
⁶³ Cu	Medium	0.6
⁶⁶ Zn	Medium	1.5
⁷⁹ Br	Medium	5.0
⁷⁵ As	High	2.0
¹⁵⁵ Gd	High	0.2

Detection Limits in 50 mg/L Uranium



²²⁶Ra, External Calibration, 1, 4, 6 pg/L, R² = 1.000




Standard Addition Analysis of Arsenic in HCl (10 % v/v)


Signal Stability and Detection Limits

Sector field ICP-MS offers superior ion transmission stability due to the high acceleration voltage (-8 kV) and excellent focusing properties.

In combination with the high sensitivity, low background, interference-free measurements and an advanced sample introduction system, the Finnigan ELEMENT2 delivers the lowest detection limits – independent of the sample matrix. The resulting low detection limits combined with highest stability at single digit ppt level enable quantification at the lowest concentrations.



1 ng/L Spike Recovery UPW



Stability Specifications

- < 1 % RSD in 10 minutes
- < 2 % RSD in 1 hour

Mass Stability

The mature design of the Finnigan ELEMENT2 delivers the fastest scan speed ever for a magnetic sector field ICP-MS. With the Finnigan ELEMENT2 a jump from 7 to 240 to 7 is realized in < 150 ms.

Scanning

Routine scanning is performed by a combination of magnetic and electric jumps.

The Finnigan ELEMENT2 has the unique capability to scan +30% from the magnet mass by decreasing the acceleration voltage. The magnetic field is kept constant while the acceleration voltage is varied. This ingenious combination of scan techniques delivers the fastest scanning ever realized by a sector field ICP-MS.

Mass Stability

The Finnigan ELEMENT2 guarantees the highest mass stability of any ICP-MS. This enables fast peak top jumping analyses rather than scanning across the whole peak, thus significantly decreasing analysis times. A combination of mature hardware and intelligent software guarantees mass stability in high resolution, which makes mass calibration a rare event.


The synthesis of these characteristics combined with high sensitivity, allowing short integration times, makes the Finnigan ELEMENT2 the fastest sector field ICP-MS ever.

The scan speed of the Finnigan ELEMENT2 in combination with with the high sensitivity opens the dimension of its use as a detector for fast transient signals.


Specifications

Mass stability: 25 ppm / 8 hours

Magnetic scan speed: m/z 7 to 240 to 7 < 150 ms




Timing for a Full Magnet Cycle




12-Hour Mass Stability for ⁷⁸Se in high resolution


Ion Detection System



Finnigan ELEMENT2 Ion Detection System



Simultaneous Measurement of the Counting and Analog Signals



The goal of real life analysis is the simultaneous determination of major and trace elements.

The Finnigan ELEMENT2 is equipped with a discrete dynode detection system that enables the quantification of both trace and major elements in a single analysis across a dynamic range of 10^9 .

The Finnigan ELEMENT2 detector system incorporates a conversion dynode at -8 kV. The high acceleration to -8 kV results in a mass independent detector response, enabling fully automatic cross calibration between the counting and analog modes.

No user interaction is necessary.

The automatic cross calibration ensures a constantly updated detector response at all times – ready for the unexpected sample!

Additionally the detector enables the measurement of large isotope ratios. With the minor isotope in counting mode and the abundant isotope in analog mode the highest precision and accuracy can be obtained.

Specifications

Dynamic Range $> 10^9$ with automatic gain calibration

The Software Suite

The Finnigan ELEMENT2 Software

The Finnigan ELEMENT2 software controls and monitors all instrument functions for ICP-MS analysis. This includes data acquisition and auto-tuning of the Finnigan ELEMENT2.

The software also provides the full range of quantification procedures required in elemental analysis (qualitative, quantitative, semi-quantitative, isotope dilution) as well as isotope ratio and time resolved analysis modes.

Controlling and Tuning

- Autotuning of all parameters, including ICP parameters, torch position, lenses and multiplier voltage
- Fully automated and configurable plasma start and stop sequence
- Easy autosampler setup using a graphical display of the autosampler

Setup of Methods

- Intuitive and easy selection of target isotopes in periodic table or spreadsheet display modes
- Automatic, customized isobaric interference correction
- Spreadsheet style 'click-and-drag' cell fill down



Creating and Running Sequences

- Intuitive and easy creation of sample analysis sequences in graphical or spreadsheet display modes
- Integrated powerful QA/QC package that meets internationally regulated requirements including US EPA 200.8, 6020. The flexible editor included can be used to define specific QA/QC criteria in any laboratory.

Displaying Results and Creation of Reports

- Real time display of spectra, calibration curves, fully quantitative results and time resolved analyses
- On-line export of time resolved data in several formats (ASCII, GRAMS, Spectacle, GLITTER, ANDI Xcalibur® for further analysis in third party programs)

- Remote Control & Diagnostic
- Setting & Read-Back of All Instrument Parameters



Instrument Startup and Tune

The Finnigan ELEMENT2 software is a state-of-the-art, simple to use software suite taking advantage of the reliability and stability of the modern operating system Microsoft® Windows® XP Professional. The software package is optimized for the needs of the routine analyst, providing stability and ease of use for basic operation of the Finnigan ELEMENT2, and yet retains the flexibility for advanced operation.


Due to the use of the Microsoft Windows XP Professional operating system and standardized programming, the data system computer can be easily connected to a network, enabling data transfer and allowing remote control of the Finnigan ELEMENT2.

Product Specification

Specifications and Installation Requirements

Sensitivity (Concentric Nebulizer)	> 1 x 10 ⁹ counts per second (cps)/ppm In	
Detection Power	< 1 ppq for non-interfered nuclides	
Dark Noise	< 0.2 cps	
Dynamic Range	> 10 ⁹ linear with automatic gain calibration	
Mass Resolution	300, 4000, 10,000 (10 % valley, equivalent to 5 % height); 600, 8000, 20,000 (FWHM)	
Signal Stability	< 1 % RSD over 10 minutes < 2 % RSD over 1 hour	
Scan Speed (magnetic)	m/z 7 to 240 to 7 < 150 ms	
Scan Speed (electric)	1 ms/jump, independent of mass range	
Oxide and Doubly Charged Ions	ratio	measured
	BaO ⁺ /Ba ⁺	< 0.002
	Ba ²⁺ /Ba ⁺	< 0.03
Power	3-phase, 230/400 V ±10 %, 50/60 Hz fused 32 A per phase Power consumption: ~9 kVA	
Environment	Temperature 18 – 24 °C (64 – 75 °F) Humidity 50 – 60 %, noncondensing, non-corrosive	
Cooling Water	~ 200 l/h Temperature 10 – 20 °C 4 – 6 bar (43 – 65 psi)	
Argon	Purity 99.996 min. 18 L/min Regulated pressure 8 –10 bar (116 –145 psi) Uninterrupted argon supply recommended	
Plasma Exhaust	1 x 6 cm Ø; 90 m ³ /h (Argon + suspended sample)	
Electronics Exhaust	2 x 15 cm Ø; 800 m ³ /h	

Finnigan ELEMENT2, Footprint and Dimensions (all dimensions in cm)



Plasma Capabilities from Thermo Electron

The use of an Inductively Coupled Plasma source (ICP) is the accepted and most powerful technique for the analysis and quantification of trace elements in both solid and liquid samples. Its applications range from routine environmental analyses to the materials industry, geological applications to clinical research and from the food industry to the semiconductor industry.

Thermo Electron Corporation is the only instrument manufacturer to offer the full range of Inductively Coupled Plasma Spectrometers (ICP-OES, Quadrupole and Sector ICP-MS) to satisfy every aspect of plasma spectrometry from routine to highly demanding research applications.

Develop your lab from the easy-to-use IRIS Intrepid II ICP-OES to the high performance X Series Quadrupole ICP-MS and up to the ultra-sophisticated Finnigan™ ELEMENT2 and Finnigan™ NEPTUNE Sector ICP-MS instruments. Each instrument combines leading-edge technology, fit for purpose and affordability with a tradition of quality, longevity, accuracy and ease of use.

IRIS Intrepid II



X Series ICP-MS



Finnigan ELEMENT2



Finnigan NEPTUNE



In addition to these offices, Thermo Electron Corporation maintains a network of representative organizations throughout the world.

Australia

+61 2 9898 1244 • analyze.au@thermo.com

Austria

+43 1 333 50340 • analyze.at@thermo.com

Belgium

+32 2 482 30 30 • analyze.be@thermo.com

Canada

+1 800 532 4752 • analyze.ca@thermo.com

China

+86 10 5850 3588 • analyze.cn@thermo.com

France

+33 1 60 92 48 00 • analyze.fr@thermo.com

Germany

+49 6103 4080 • analyze.de@thermo.com

Italy

+39 02 950 591 • analyze.it@thermo.com

Japan

+81 45 453 9100 • analyze.jp@thermo.com

Netherlands

+31 76 587 98 88 • analyze.nl@thermo.com

Nordic

+46 8 556 468 00 • analyze.se@thermo.com

South Africa

+27 11 570 1840 • analyze.sa@thermo.com

Spain

+34 91 657 4930 • analyze.es@thermo.com

Switzerland

+41 61 48784 00 • analyze.ch@thermo.com

UK

+44 1442 233555 • analyze.uk@thermo.com

USA

+1 800 532 4752 • analyze.us@thermo.com

www.thermo.com



Thermo Electron (Bremen) GmbH is certified DIN EN ISO 9001:2000

©2004 Thermo Electron Corporation. All rights reserved. Microsoft and Windows are registered trademarks of the Microsoft Corporation. All other trademarks are the property of Thermo Electron Corporation and its subsidiaries.

Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

BR30027_E 05/04C

Thermo
ELECTRON CORPORATION