# Thermo Scientific Nicolet iS50 FT-IR Spectrometer

The Materials Analysis Workstation

The Thermo Scientific Nicolet iS50 FT-IR spectrometer uniquely combines multi-tasking capabilities and high performance in an affordable, optimized footprint system. Extensive upgrade options include a built-in ATR, an automated beamsplitter exchanger and Raman, NIR, TGA-IR and GC-IR modules.





The main goal of the busy analytical laboratory, whether working with polymers, rubbers, pharmaceuticals, forensics or any other materials, is answering specific questions. The Nicolet™ iS™50 FT-IR spectrometer provides smarter tools — both in the instrument and in the software — to lead you to definitive answers.

The Nicolet iS50 couples multiple sources, sampling stations and detection options through Touch Point one-touch setup and operation, especially when driven by the iS50 ABX automated beamsplitter exchanger. Built-in tools leave open options, such as simultaneous installation of a sample compartment iS50 Raman module and the iS50 ATR multi-range, diamond sampling station.

All offerings you have come to know, like Thermo Scientific Smart Accessories and System Performance Verification, are supported on the Nicolet iS50.

In the most demanding laboratory situations, full validation is available along with hyphenated tools like the iS50 GC-IR module and the TGA-IR accessory. The Thermo Scientific OMNIC software contains new tools such as the Mercury TGA and Mercury GC analysis routines, auto-reporting and the archiving of analysis results when you save your data.

Experience FT-IR beyond the ordinary with the Nicolet iS50 FT-IR spectrometer.



## Go beyond your expectations

with the Nicolet iS50 FT-IR Spectrometer



## **Thermo Scientific** Nicolet iS50 FT-IR Spectrometer

The Materials Analysis Workstation

#### **Beam Size and Optical Filters**

Continuously variable J-stop

**High Efficiency Sealing System**  Sealed and desiccated standard Purge connections standard

> **Motorized ZnSe Wire** Grid Polarizer, in/out and

rotation controlled

Small volume

- High-resolution standard, 0.09 cm<sup>-1</sup>
- Filter wheel accepts up to five industry standard one inch filters for visible, far-IR etc.

#### **Four Position Source Mirror**

- Polaris<sup>™</sup> Long-lifetime mid-IR source
- Tungsten-Halogen NIR/Vis source
- · Raman InGaAs detector
- Focused emission port

Source	High (cm <sup>-1</sup> )	<b>Low</b> (cm <sup>-1</sup> )
Polaris IR	9600	20
NIR/Vis	27,000	2000
External	Custom	Custom

### **Three Position Detector Mirror**

- User replaceable, LN2 cooled
- DLaTGS (standard)
- User replaceable, room temperature

Detector	High (cm <sup>-1</sup> )	<b>Low</b> (cm <sup>-1</sup> )
DLaTGS-KBr	12,500	350
MCT-High D*	11,700	800
MCT-A	11,700	600
MCT-B	11,700	400
Time-resolved MCT	11,700	650
Silicon	27,000	8600
PbSe	11,000	2000
InGaAs	12,000	3800
InSb	11,500	1850
DLaTGS-Csl	6400	200
DLaTGS-Polyethylene	700	50
Si bolometer	600	15
Photoacoustic	10,000	400

### **Full Sized Sample Compartment**

- KBr or Csl windows
- Motorized purge shutters
- Compatible with standard and Smart Accessories

#### **Small Footprint**

System	Weight	<b>Dimensions</b> (W $\times$ D $\times$ H)
Base	60 kg 132 lbs	$62.6 \times 69.8 \times 27.6$ cm $25 \times 27 \times 11$ in
With ABX	64 kg 141 lbs	$62.6 \times 69.8 \times 50.8 \text{ cm}$ $25 \times 27 \times 20 \text{ in}$

- No added footprint for Raman and dedicated ATR
- Only 27.9 cm (11 in) more with iS50 NIR module

#### **Optional iS50 ABX Automated Beamsplitter Exchanger**

- Up to three beamsplitters
- Cover far-mid-near, far-mid-vis or other combinations
- Less than 25 seconds per exchange
- 52 mm beamsplitter diameter



#### **Multiple Beamsplitter Options**

Beamsplitter	High (cm <sup>-1</sup> )	Low (cm-1)
KBr	7800	350
XT-KBr	11,000	375
Quartz	27,000	2800
CaF <sub>2</sub>	14,500	1200
Csl	6400	200
Solid-Substrate	700	20

Csl is offered as a dedicated system

#### **Sample Compartment**

#### Fluorescence-Free FT-Raman

- 1064 nm diode laser
- Weight: 7.6 kg (16.8 lbs)
- Full validation for regulated environments available
- Built in x-y-z stage
- Point and collect, area map, well plates
- Driven by μView, Atlμs and Array Automation
- Screening and cluster analysis



#### **Dynamically Aligned Interferometer**

- Thousands of field-proven installs
- Durability and speed
- Tilt and shear full mirror control

#### **Easy Laser Replacement**

- Modular design
- Externally mounted

#### **Validation / Attenuation Wheel Standard**

- NIST traceable 1.5 mil polystyrene
- NIST traceable NG-11 glass

Far-infrared spectrum of acetylferrocene

• Two selectable energy screens



#### **Optional iS50 ATR**

- Built-in, all-reflective diamond ATR
- Mid- to far-IR capable: 80 to over 5000 cm<sup>-1</sup>
- · Monolithic diamond for durability
- Software-controlled activation
- Pressure applied to 60 lbs
- · Removable tray for cleaning
- Liquid/volatiles cover available
- Full validation for regulated environments available



# iS50 Raman Module **Built-in USB Video Camera** • 5 mm field of view View and collect Archive images **Quick Lock Connections** Cable free installation **External iS50 NIR Module**

#### **Integrating Sphere and NIR Fiber Port**

- Supports sample cup spinner, viscous sample accessory, and more
- Fiber optic SMA in/out ports
- Optional Thermo Scientific SabIR probe
- Full validation for regulated environments
- Available in integrating sphere only configuration

#### **External iS50 GC-IR Module**

#### **Heated Flow Cell and Transfer Line**

Raman map and video image of fiber

• 300 °C continuous operation



Infrared data management and autosampler control when equipped with Thermo Scientific TRACE 1310 GC and Thermo Scientific Dionex Chromeleon chromatography data system software.

#### **Nicolet iS50 Specifications**

Spectrometer	
Polaris High Stability, Long Lifetime Mid-IR Source	Standard
Tungsten-Halogen Near-IR/Visible Source	Option
Four Position Source Mirror	Option
Continuously Variable Iris Aperture	Standard
Gold Optical Coatings	Standard
Aluminum Optical Coatings	Option
DLaTGS Detector	Standard
Three Position Detector Mirror	Option
Attenuation Wheel	Standard
Validation Wheel	Standard
Automated Polarizer	Option
Automated Filter Wheel	Option
Automated Beamsplitter Exchanger	Option
Automated Sample Compartment Purge Shutters	Option
A/D Converter	24 bit
Interface	USB 2.0
Software	
Operating System	Windows® 7
OMNIC Software	Standard
ValPro System Validation Software	Option
21 CFR Part 11 Compliance Tools	Option
External Beam Capabilities	
Dual Side External Beams	Option
Collimated Emission Port	Option
Focused Emission Port	Option
Side External Detector Port	Option

Performance Specifications	
Spectral Range, Standard System	7800-350 cm <sup>-1</sup>
Spectral Range, Csl Optics	6400-200 cm <sup>-1</sup>
Spectral Range, Multi-Range Optics	27,000-20 cm <sup>-1</sup>
Optical Resolution, Mid-IR	Less than 0.09 cm <sup>-1</sup>
Signal-to-Noise, 1 minute scan, Peak-to-Peak, 4 cm <sup>-1</sup>	55,000:1
Signal-to-Noise, 5 second scan, Peak-to-Peak, 4 cm <sup>-1</sup>	13,000:1
Ordinate Linearity	0.07%T
Wavenumber Precision	Better than 0.01 cm <sup>-1</sup>
Scan Velocity (15 values)	0.158-6.28 cm/sec
Rapid Scan, Spectra Per Second	65 (at 16 cm <sup>-1</sup> ), 95 (at 32 cm <sup>-1</sup> )
MCT Dewar LN₂ Hold Time	18 Hours
Physical Characteristics	
Spectrometer Weight	60 kg (132 lbs)
Spectrometer Dimensions $(W \times D \times H)$	$62.6 \times 69.8 \times 27.6 \text{ cm}$ $25 \times 27 \times 11 \text{ in}$
Sample Compartment Dimensions $(W \times D \times H)$	$21 \times 26 \times 15 \text{ cm}$ $8.3 \times 10.2 \times 5.9 \text{ in}$

5 Years 1 Year

**(**E **(**B)

Performance Specifications

## Optional Application Modules

In addition to the application modules listed below, a full line of accessories are available, including: infrared microscopes, TGA-IR, and Linear-scan Dual-channel Collection. The Nicolet iS50 is compatible with both standard and Smart Accessories.

iS50 ATR Module	
Crystal	Diamond
Spectral Range	5000-80 cm <sup>-1</sup>
Down Force of Pressure Device	60 lbs
Detector	Proprietary DLaTGS
iS50 Raman Module	
Laser	1064 nm diode
Laser Power	>450 mW at sample
Laser Spot Size	60 microns
Sampling Plates	48 well, 9 well, vials, microscope slides
Stage Resolution	5 micron steps
Weight	7.6 kg (16.8 lbs)
Compliance	Class 1 laser product

ntegrating Sphere Window	Sapphire
nterior of Integrating Sphere	Gold coated
ntegrating Sphere Detector	InGaAs
Fiber Optic Connections	Standard SMA

Fiber Optic Connections Standard SM
Fiber Optic Detector InGaAs
Validation Wheel Standard

Mid-infrared Source and Interferometer Warranty

## Gas Cell

iS50 NIR Module

**Other** 

Spectrometer Warranty
Regulatory Approvals

Gas Cell	15 cm $\times$ 1 mm gold-coated light pipe
Temperature	300 °C max transfer line and cell heaters, USB controlled
Detector	LN <sub>2</sub> -cooled MCT-A
Exhaust Line	Passes through activated charcoal filter to rear panel fitting

#### thermoscientific.com

©2012 Thermo Fisher Scientific Inc. All rights reserved. Windows is a registered trademark of Microsoft Corporation. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Africa +27 11 822 4120 Australia +61 3 9757 4300 Austria +43 1 333 50 34 0 Belgium +32 53 73 42 41 Canada +1 800 530 8447 China +86 10 8419 3588

Denmark +45 70 23 62 60 Europe-Other +43 1 333 50 34 0 Finland/Norway/Sweden +46 8 556 468 00 France +33 1 60 92 48 00 Germany +49 6103 408 1014

India +91 22 6742 9434 Italy +39 02 950 591 Japan +81 45 453 9100 Latin America +1 561 688 8700 Middle East +43 1 333 50 34 0 Netherlands +31 76 579 55 55 Thermo Electron Scientific Instruments LLC, Madison, WI USA is ISO Certified. Class I Laser product.

New Zealand +64 9 980 6700 Russia/CIS +43 1 333 50 34 0 Spain +34 914 845 965 Switzerland +41 61 716 77 00 UK +44 1442 233555

USA +1 800 532 4752



Part of Thermo Fisher Scientific