

Comparing IC and HPLC systems for speciation analysis – a case study

Authors: Daniel Kutscher¹, Carl A. Fisher², and John Schmelzel³; ¹Thermo Fisher Scientific, Bremen, Germany; ²Thermo Fisher Scientific, Sunnyvale, CA, USA; ³Thermo Fisher Scientific, San José, CA, USA

Keywords: ion chromatography, high performance liquid chromatography, chromium (III) and (VI), IC-ICP-MS

Goal

This technical note compares the performance of ion chromatography (IC) and high performance liquid chromatography (HPLC) methods for speciation analysis of chromium, including method related benefits and drawbacks.

Introduction

Speciation analysis is a widely accepted tool for the quantitative determination of different chemical forms, or species, of an element. Two of the primary examples of the need to perform speciation analysis are chromium¹ and arsenic^{2,3}, as both elements may appear as different species (e.g., Cr (III) and Cr (VI) or inorganic and organic arsenic), each with different toxicity and thus, potential related hazards. Another element that requires speciation is bromine, as ozonation of drinking water may lead to the formation of carcinogenic bromate (BrO_3^-) as a byproduct⁴.



The different species of arsenic are typically charged depending on the pH, therefore ion exchange chromatography is the preferred technique for separation. In most cases, anion exchange chromatography using an ammonium carbonate-based eluent allows complete separation with narrow peak shapes and excellent signal-to-noise ratios. For chromium, the two relevant species are charged differently, so that, ideally, a mixed-mode resin is used containing separation capacity for both cationic and anionic species. The separation of bromide and bromate (among other brominated organic molecules) can be accomplished using anion exchange chromatography, using, for example, KOH as an eluent (combined with post-column eluent suppression).

IC systems are designed to separate charged molecules, and are, therefore, the separation instruments of choice for all of the elements mentioned above. As an additional benefit for ICP-MS detection, background contamination is minimized as IC systems have a metal-free flow path.

For Cr speciation, however, other methods, such as ion pairing chromatography, can be applied. As this chromatographic method depends on the use of complexing agents, such as EDTA and organic modifiers, significant spectral interferences will be created affecting the main isotope of chromium, ^{52}Cr . Typical elution conditions for IC and HPLC-based methods for separation of both chromium species of interest are summarized in Table 1. The table also contains information about potential polyatomic interferences of the commonly used isotopes of chromium, which lead to an increased background signal in ICP-MS. Moreover, not all compounds required to prepare the eluent may be available in sufficient purity.

Table 1. Eluent composition, elution conditions, and resulting interferences for IC- and HPLC-based methods for speciation analysis of chromium

Method	Eluent	Elution conditions	Potential interferences
IC (Ion exchange)	2% NH_4NO_3	Isocratic	None
HPLC (Ion pairing chromatography)	EDTA in combination with different buffers, such as TBAP in methanol/water	Isocratic or gradient	$^{40}\text{Ar}^{12}\text{C}^+$ for $^{52}\text{Cr}^+$ $^{40}\text{Ar}^{13}\text{C}^+$ for $^{53}\text{Cr}^+$

The choice of the right separation method for a given analytical problem often determines the most appropriate hardware, such as type of chromatography system, column dimensions, stationary and mobile phase. Whereas most applications are based on liquid chromatography (Gas chromatography (GC)-ICP-MS is an alternative for volatile species or species that have volatile derivatives⁵), there is a significant difference between IC and HPLC methods, even though the required hardware can seem almost identical.

Ion chromatography

- Hardware does not use parts made of metal-containing components (e.g., stainless steel or titanium alloys) that come into contact with the eluent. This results in significantly reduced backgrounds due to the absence of metal leaching over time, and assures compatibility with strong eluents such as diluted acids or alkaline eluents, such as hydroxides or carbonates.
- To assure best compatibility with the ICP-MS sample inlet, eluents such as KOH can be neutralized to water using a suppressor before introduction to the plasma source.
- Modern IC systems can deliver elevated pressures (up to 6000 psi) to permit use of higher resolution columns that yield excellent peak separation and short elution time.

High performance liquid chromatography

- Hardware often uses tubing or connectors made of stainless steel or pump heads made of titanium. Bio-inert systems contain significantly fewer parts made of metal to reduce unwanted effects, such as unfolding or degradation of proteins, but do not achieve the low backgrounds observed with IC.
- Some separations will require the addition of ion pairing agents in order to achieve charge separation or neutralization (example: Cr speciation). This will lead to increased backgrounds through interferences.

Experimental

Materials and methods

A Thermo Scientific™ iCAP™ RQ ICP-MS system was used for determination of background levels in different eluents after they had passed through either an IC or an HPLC system. The instrument was tuned on each measurement day to achieve the best sensitivity across the mass range and effective suppression of occurring spectral interferences. For IC, a Thermo Scientific™ Dionex™ Integriion™ HPIC™ system was used. The IC system was compared against a standard HPLC system. Please note that the use of this eluent may be outside of working specifications of typical HPLC systems (typically between pH 1-13). The time-dependent evaluation of eluent contamination was accomplished either online or offline (fraction collection).

It is possible that this contamination may have resulted from the process of fraction collection or from the sample containers used in this study. Further investigation will be directed to determine the origin of the observed background. For some elements, a time-dependent behavior of elution was observed for the HPLC system under test in this study, leading to a reduction in the observed levels over time. Generally, metal alloys are known for forming a passive oxide layer on the surface that protects the core of the metal against corrosion when in contact with certain solutions and chemicals. For stainless steel, the passive oxide layer mainly consists of a Cr-oxide, but also Fe and Ni can be present in different ratios. The ratio and the depth of the oxide layer depends on the properties of the solution (e.g., pH, oxidizing/reducing potential, temperature, duration) with which it is in contact. This indicates that the passivating layer has not reached a steady state, or in other words, the HPLC system has not yet been fully equilibrated with the eluent in the 150 min of this study. The dependency of the passivating layer from the composition of the eluent, however, means that not only a leaching of metals is likely to be observed when running an IC method on an HPLC system, but also a potential change of the leaching behavior can be expected when switching to a different eluent/mobile phase for other applications. However, in no case did the observed levels of these elements reach the levels observed with a dedicated IC system. As an example, the time-dependent elution of chromium is shown in Figure 1.




Figure 1. Time-dependent elution of ⁵²Cr as a contaminant in an IC and HPLC system

Conclusion

Eluent contamination with trace elements is significantly increased for a variety of elements on an HPLC system compared to an IC system. This is most likely caused by metal-containing parts in the flow path of the eluent/mobile phase, including the pump heads, connectors, and tubing. For some elements such as chromium or vanadium, a time-dependent wash out was observed in the first two hours of operation; however, equivalent background levels between the IC and HPLC systems were not achieved. Systemic accumulation of metal contamination after long-term usage of an HPLC system with a typical IC eluent was not evaluated here, but, based on the current findings, it would likely be significant, ultimately degrading system performance. To avoid elevated backgrounds caused by leaching of metals from hardware components, chromatographic methods based on ion-exchange chromatography should be run on a dedicated IC system.

References

1. Thermo Fisher Scientific, Application Note 44407: Determination of chromium species using ion chromatography coupled to inductively coupled plasma mass spectrometry. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/an-44407-ic-icp-ms-chromium-drinking-water-an44407-en.pdf>
2. Thermo Fisher Scientific, Technical Note 43357: Arsenic Speciation in human urine by hyphenated ion chromatography (IC) and inductively coupled plasma mass spectrometry (ICP-MS). <https://assets.thermofisher.com/TFS-Assets/CMD/Technical-Notes/TN-43357-IC-ICP-MS-Arsenic-Urine-TN43357-EN.pdf>
3. Thermo Fisher Scientific, Application Note 43255: Determination of inorganic arsenic in rice using IC-ICP-MS. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-43255-IC-ICP-MS-Inorganic-Arsenic-Rice-AN43255-EN.pdf>
4. Thermo Fisher Scientific, Application Note 43227: Speciation of bromine compounds in ozonated drinking water using ion chromatography and inductively coupled plasma mass spectrometry. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-43227-IC-Bromine-Water-AN43227-EN.pdf>
5. Thermo Fisher Scientific, Application Note 43371: Ultra trace tin speciation with GC-ICP-MS using the GCI-100 interface. <https://assets.thermofisher.com/TFS-Assets/CMD/Application-Notes/AN-43371-GC-ICP-MS-Ultratrace-Tin-Speciation-AN43371-EN.pdf>
6. Specification Sheet Thermo Scientific™ Dionex™ IonPac™ AS7 IC Column; Thermo Fisher Scientific

Find out more at thermofisher.com/speciationanalyzer

ThermoFisher
SCIENTIFIC