

Discover Bio

Microwave Peptide Synthesis

CEM's best-selling line of microwave peptide synthesizers now includes the new Discover Bio System, a fully-featured manual microwave peptide synthesizer that offers the most affordable option for adding microwave technology to your peptide synthesis laboratory.

Discover Bio is a compact, single-mode microwave system with a liquid handling module for the transfer of amino acids, reagents and solvents. Temperature is controlled using a fiber optic probe for the most accurate measurement.



Discover Bio System

Manual Microwave Peptide Synthesizer

- Greater purity peptides than conventional synthesis
- Make peptides faster than ever with cycle times under 15 minutes
- Replaceable frit provides low consumable cost
- Flexible system can be used for both peptide & small molecule organic synthesis
- Upgradable to fully automated Liberty Blue system
- License included allowing the use of microwave energy for both the coupling & deprotection steps

CEM
We are peptide synthesis



- Integrated washing and reagent delivery system
- Two reaction vessels cover wide scale range (0.05 to 1 mmol)
- Fiber optic probe for the most accurate temperature measurement
- Inert atmosphere provides cleaner chemistry
- Nitrogen bubbling for effective agitation
- Vessel's integrated spray head allows efficient washing

Best-in-Class Applications Assistance & Technical Support

Should you have need of service or applications support, CEM maintains a fully-trained staff of chemists and service technicians to ensure that your questions are answered in a timely manner. We offer technical service support 24/7 and our peptide chemists are just a phone call or e-mail away. That's why prestigious research facilities, universities, and major pharmaceutical and biotech companies around the world turn to CEM for their peptide synthesis needs.

CEM Corporation

PO Box 200
 Matthews, NC 28104
 3100 Smith Farm Road
 Matthews, NC 28106
 Tel: 800-726-3331
 Tel: 704-821-7015
 Fax: 704-821-7894
 Email: info@cem.com
www.cem.com

France

CEM μ Wave S.A.S.
 Immeuble Ariane
 Domaine Technologique de Saclay
 4, rue René Razel
 91892 ORSAY Cedex
 Tel: +33 (01) 69 35 57 80
 Fax: +33 (01) 60 19 64 91
 Email: info.fr@cem.com
www.cemfrance.fr

Germany

CEM GmbH
 Carl-Friedrich-Gauss-Str.9
 47475 Kamp-Lintfort
 Tel: (49) 2842-9644-0
 Fax: (49) 2842-9644-11
 Email: info@cem.de
www.cem.de

Ireland

CEM Technology (Ireland) Ltd.
 Sky Business Centre
 9a Plato Business Park
 Damastown
 Dublin 15
 Tel: +353 (0) 1 885 1752
 Fax: +353 (0) 1 885 1601
 Email: info.ireland@cem.com
www.cemmicrowave.co.uk

Italy

CEM S.R.L.
 Via Dell' Artigianato, 6/8
 24055 Cologno al Serio (Bg)
 Tel: (39) 35-896224
 Fax: (39) 35-891661
 Email: info.srl@cem.com
www.cemmicroonde.it

Japan

CEM Japan K.K.
 2-18-10 Takanawa
 Minato-ku, Tokyo
 108-0074
 Tel: +81-3-5793-8542
 Fax: +81-3-5793-8543
 Email: info@cemjapan.jp
www.cemjapan.co.jp

United Kingdom

CEM Microwave Technology Ltd.
 2 Middle Slade
 Buckingham Industrial Estate
 Buckingham MK181WA
 Tel: (44) 1280-822873
 Fax: (44) 1280-822342
 Email: info.uk@cem.com
www.cemmicrowave.co.uk