

NEW



Disposable Teflon Liners for MARSXpress Vessels

Save time and money.

The Teflon® liner protects the vessel from batch-to-batch contamination without washing vessels between runs. Simply insert a disposable liner into the MARSXpress™ vessel, cap it, and run. When the digestion is complete, transfer your solution to an appropriate flask or vial and dispose of the liner. You are ready to start your next sample. Disposable liners are the perfect accessory for high-throughput labs running USEPA 3051A and 3015A methodology. They can also be used for any digestion application run at or below 180 °C for 30 minutes or less.



We Simplify Science

cem.com

Get clean digests and eliminate vessel cleaning.

Disposable Teflon liners are clean enough to achieve < 5 ppb recoveries when run with acid only.

Element	Concentration (ppb)
Ca	8
Al	<5
Na	12
Ba	<5
Cu	<5
Fe	<5
Pb	<5
Mn	<5
Ni	<5
Bi	<5
Co	<5
V	<5
Cr	<5
Sr	<5
Se	<5
Zn	<5
As	<5
Cd	<5



For more information or to order, contact CEM Corporation at 704-821-7015 or info@cem.com.

Liner Forming Kit

Part: 404450

Includes: 1 Flaring Tool
and 1 Insertion Tool

20 mL Disposable Liners

Part: 404495

Quantity: 100 Liners

55 mL Disposable Liners

Part: 404460

Quantity: 100 Liners

75 mL Disposable Liners

Part: 404465

Quantity: 100 Liners

Note: iWave® Temperature Sensing required. Contact CEM to upgrade.